
Microarray Data Analysis (1)

Data quality
assessment
for microarrays

Microarray Raw Data (1): cDNA arrays

One .GPR (for GenePix Results) file per chip containing

- One row per gene but many columns with
 - Intensity values for each channel (R, G)
 - Summary values for intensities
 - Quality controls, such as FLAG

Intensity values are converted into a single expression matrix containing

- One column per chip with $\log(R/G)$ values
- One row per gene (same rows as .GRP files)
- Gene information stored in a .GAL (for GenePix Array List) file
- Both .GPR and .GAL are ASCII files
- An accurate description of these files is available [here](#)

Microarray Raw Data (2): Affymetrix arrays

One .CEL file per chip, containing

- PM and MM values for each probe in the chip
- Presence/Absence calls (one per probeset)

They can be interpreted as a statistical test of the spot foreground intensity in the experimental sample respect to the background intensity distribution

- Separate PM/MM values are converted into a single expression matrix containing:
 - One column per chip with absolute intensity values
 - One row per probeset
- Gene information stored in the .CDF file
- .CEL file is a binary file
- An accurate description of these files is available [here](#)

Looking at microarray data Diagnostic Plots

*Was the experiment a
success?*

Exploring experimental results

- Microarray experiments generate huge quantities of data.
- It is hard to decide if things “seem to be all right” just by looking at the numbers.
- Standard statistical approach: use plots.
 - Show all data together.
 - Highlight structures,
 - May help detect problems (“unusual patterns”)

Diagnostic plots for microarrays

- Microarray data usually considered at two levels
 - Low-level: Data directly coming from the scanner
 - High-level: Processed from low-level data. Expression values, normalized or not.
- Some diagnostic plots may differ between one and two color arrays, specially for looking at low level values.
- Other may be used for any type of arrays or for any level.


Diagnostic plots

<i>Microarray type</i> <i>Data type</i>	One color	Two color	General
Low level 1col: (probe, probe-set) 2col: (single channel)	Layout image Degradation plots Density plots <i>Probeset plots</i>	Scatterplots R,G MA-plot Signal2Noise plots Layout image (G, R)	PCA Histogram/Density Boxplot
High-level 1 col: Relative expression 2 col: Absolute expression	MA-Plots <i>Model-based plots</i> (NUSE, RLE, Residual)	Layout image (log ratios)	PCA Histogram/Density Boxplot


Diagnostic plots for two color arrays

Red / Green overlay images

- Start by looking at the slides



Bad: high bg



Good: low bg

Provides information on

- colour balance,
- uniformity of hybridization,
- spot uniformity,
- background, and
- artifacts such as
 - dust or
 - scratches


MA-plot (1)

To determine whether correction (normalization) is needed, one can plot R vs G intensities and see whether the slope of the line is around 1


But a better representation of genes with “medium” expression is to take logs...

MA-plot (2)

Biologically, a unit change in \log_2 represents a 2-fold change



Linear



Log scale

MA-plot (3)

- An improved method of the R vs G plot is the *MA-plot*, which is basically a scaled 45 degree rotation
- It is a plot of the distribution of
- **M-value** which is the \log_2 of the R/G intensity ratio


$$M = \log_2(R/G) = \log_2(R) - \log_2(G)$$

A-value which is the \log_2 of average intensity


$$A = (\log_2(R) + \log_2(G))/2 = \log_2 \sqrt{RG}$$

- The general assumption is that most of the genes would not see any change in their expression -> the majority of the points on the M would be located at 0, since $\log_2(1) = 0$


MA-plot (4): M vs A



MA-plot (5)





$\log_2 R$ vs $\log_2 G$









$M(=\log_2 R/G)$ vs $A(=\log_2 \sqrt{RG})$

MA-plot for spotted arrays (2 colors)




Signal/Noise histograms



Images with high background tend to have lower $\log_2(\text{signal/noise})$ ratios



Quality between slides





Diagnostic plots for affy chips

Image plots for affymetrix chips


D
O
O
G




B
A
D



MA-plot for GeneChip arrays (1 color)



M
 $\text{Log}_2(\text{MT}) - \text{Log}_2(\text{WT})$





A
 $\text{Log}_2(\text{MT} * \text{WT}) / 2$
(signal strength)

Density plots (1)

- Density plots of probe intensities are useful to visualize differences in distribution arrays
- Some hints for the interpretation:
 - Density is skewed to right if genes have high expression values
 - If the shape of the bell has a sharp central point and fat tails there are genes with high expression values
 - If the shape of the bell has a flat central point and skinny tails there are few genes with high expression values



Density plots (2)





Digestion or (degradation) plots (1)

- RNA degradation plots show the mRNA average expression from 5' to 3'
- They allow to assess mRNA quality from biological samples that has been used to perform arrays
- Each curve represents a single array
- Ideally, curves should be flattened as much as possible

Digestion or (degradation) plots (2)




Box plots



This plot can be used in both one and two-color arrays

Principal components (PCA) plot



This plot can be used in both one and two-color arrays